

*****SAMPLE*****

Policy and Procedure: Appropriate handling and transport of contaminated equipment and/or instruments

Purpose: To prepare contaminated instruments for transport:

1. All individuals with potential contact with the instruments have been given the required Blood Borne Pathogens training.
2. The instruments are sprayed with an enzymatic foam. Point-of-use cleaning is required to remove blood/body fluids prior to transport.
3. The contaminated instruments are placed into a puncture resistant, plastic container with a lid that is marked with an approved biohazard label.
4. The container is then transported to _____ for processing.